

Customer Information

14.04.2022, Version 1, FRI


Importing the tube test COD 60 in saltwater to the *NANOCOLOR*[®] VIS II and *UV*/_{VIS} II

Overview	
Photometer	<i>NANOCOLOR</i> [®] VIS II and <i>UV</i> / _{VIS} II
Availability	14.04.2022
Test number	0-20
Available on	Download method file
Required accessories	USB drive

Important note

Although importing of the method files does not affect the data and the costumers special methods on the instrument, MACHEREY-NAGEL recommends precautions to avoid data loss. Perform a backup of the instrument database before the import of the method file, via *Main Menu* → *Settings* → *Backup* → *USB drive* → *Backup*. Furthermore, it is recommended to save user defined methods via *Methods* → *Special Methods* → *User Methods* → *List* → *Export* to a USB drive prior to the import process. Detailed data backup instructions can be found in the manual of the respective photometer.

Test import

- Download the method ZIP-file via the link above and extract the method files on your computer.
- Copy the .mtd files to the upper level of an USB drive.
- Connect the USB drive to the photometer.
- Enter the list of the user methods via *Methods* → *Special Methods* → *User Methods* → *List*.
- Select the  -icon and import the method file for the respective photometer.
- The successful import is indicated by a message on the screen and must be confirmed with *OK*. ATTENTION: The test is not displayed in the open list view after importing!
- The *NANOCOLOR*[®] CSB 60 in saltwater tube test can be selected by entering the number 020 in the *Test number* menu or via inserting of the barcoded test tubes into the photometer.

Modifications

- Import the method file VIS_II_0_20.mtd for the implementation of the *NANOCOLOR*[®] COD 60 test in saltwater (REF 985020) on the *NANOCOLOR*[®] VIS II.
- Import the method file UVVIS_II_0_20.mtd for the implementation of the *NANOCOLOR*[®] COD 60 test in saltwater (REF 985020) on the *NANOCOLOR*[®] UV/VIS II.

Contact

For further questions, please do not hesitate to contact us:

Technical Support and Customer Service	Tel.:	+49 24 21 969-332
	E-Mail:	csc@mn-net.com
For orders (international)	E-Mail:	sales@mn-net.com