

Filter Papers – Filtrierpapiere

Grade / Sorte	MN 616 LSB - 50
Capacity / Kapazität [meq/cm²]	0,03 ± 0,005
Capacity / Kapazität [meq/Filter with ø 48 mm]	0,55 ± 0,1
Temperature stability / Temperaturstabilität	Max. 60°C
Applications and properties / Anwendungen und Eigenschaften	Filter paper with strongly basic anion exchange resin, Matrix polystyrene crosslinked with DVB, functional groups quarternary ammonium compounds, supplied in OH ⁻ form/ Filterpapier mit stark basischen Anionenaustauscherharz auf Polystyrolbasis mit DVB quervernetzt, aktive Gruppen quaternäre Ammoniumverbindungen, Lieferform OH ⁻

**These data are typical values and must not be regarded as specification.
Die angegebenen Daten sind typische Werte und stellen keine Spezifikation dar.**

