

Filter Papers – Filtrierpapiere

Grade / Sorte	MN QF-10
Weight / Gewicht [g/m²]	85
Filtration speed / Filtriergeschwindigkeit [s]	<21
Surface / Oberfläche	Smooth / glatt
Applications and properties /	Quartz microfibre filter without binder, very high purity, for air sampling /
Anwendungen und Eigenschaften	Quarzfaserfilter, binderfrei, sehr rein, für Luftanalytik /
Typical content of trace elements / typische Blindwerte [mg/kg]	Al <110, Ba <100, Ca <100, Cd <1, Co <5, Cr <10, Cu <5, Fe <50, Mg <20, Na <90, Ni <10, Pb <10, V <10, Zn <10

**These data are typical values and must not be regarded as specification.
Die angegebenen Daten sind typische Werte und stellen keine Spezifikation dar.**